

ECOTAXA

HOW TO ANALYSE THE ECOTAXA TABLES EXPORTED FROM A ZOOSCAN/ZOOPROCESS PROJECT?

Quantitative Imaging Platform of Villefranche sur Mer (PIQv)

Elineau&Picheral

Date	Name operator	Object	Chapter
2018/12/14	Elineau Amanda	creation	

Summary

1. Variables to utilize from the tables 2
 - 1.1 Variables to utilize in **Object details** (column called "object_XXXX") 2
 - 1.2 Variables to utilize in **Sample details** (column called "sample_XXXX") 2
 - 1.3 Variables to utilize in **Process details** (column called "process_XXXX") 2
 - 1.4 Variables to utilize in **Acquisition details** (column called "acq_XXXX") 2
2. Equations..... 2
 - 2.1 Concentration = Number of individus in the sampling/m³ 2
 - 2.2 Biovolume = Volume biomass of individus in the sampling/m³ 3

IMPORTANT NOTE:

Each lines of the table is 1 object (living or not living) scanned by the ZooScan
Each columns of the table is 1 variable

1. Variables to utilize from the tables

1.1 Variables to utilize in Object details (column called "object_XXXX")

object_lat = GPS latitude of the sampling (decimal format)
object_lon = GPS longitude of the sampling (decimal format)
object_annotation_status = status of the images (unclassified, dubious, predicted or validated)
object_annotation_category = taxonomic name of the images
object_area_exc = surface area of the objects excluded the empty spaces (holes) in pixel
object_area = surface area of the objects included the empty spaces (holes) in pixel
object_major = primary axis of the best fitting ellipse for the object in pixel
object_minor = secondary axis of the best fitting ellipse for the object in pixel

1.2 Variables to utilize in Sample details (column called "sample_XXXX")

sample_tot_vol = initial volume filtered = volume of seawater filtered by the net (m³)
sample_tow_nb = number of tow in the analyzed sample

1.3 Variables to utilize in Process details (column called "process_XXXX")

process_particle_pixel_size_mm

1.4 Variables to utilize in Acquisition details (column called "acq_XXXX")

acq_id = origin of the object
acq_sub_part = sample fraction analysed into the ZooScan

2. Equations

IMPORTANT NOTE: You have to convert the pixel variables in metric form

Area (mm²) = object_area x (process_particle_pixel_size_mm)²
Area excluded (mm²) = object_area_exc x (process_particle_pixel_size_mm)²
Major (mm) = object_major x process_particle_pixel_size_mm
Minor (mm) = object_minor x process_particle_pixel_size_mm

2.1 Concentration = Number of individus in the sampling/m³

Number of recurrence for each taxonomic group in the variable
object_annotation_category:

Concentration = nb. ind./m³ = (object_annotation_category x acq_sub_part) / sample_tot_vol

2.2 Biovolume = Volume biomass of individus in the sampling/m³

2.2.1 Plain biovolume

Radius of a circle = r (mm) = $\sqrt{\text{Area (mm}^2\text{)} / \pi}$

Spherical Volume = V (mm³) = $\frac{4}{3} \times \pi \times r^3$

Biovolume = Bv (mm³/m³) = (Spherical Volume x acq_sub_part) / sample_tot_vol

2.2.2 Riddled biovolume

Radius of a circle = r (mm) = $\sqrt{\text{Area excluded (mm}^2\text{)} / \pi}$

Spherical Volume = V (mm³) = $\frac{4}{3} \times \pi \times r^3$

Biovolume = Bv (mm³/m³) = (Spherical Volume x acq_sub_part) / sample_tot_vol

2.2.3 Ellipsoide biovolume

Spherical Volume = V (mm³) = $\frac{4}{3} \times \pi \times [(\text{Major(mm)}/2) \times (\text{Minor(mm)}/2) \times (\text{Minor(mm)}/2)]$

Biovolume = Bv (mm³/m³) = (Spherical Volume x acq_sub_part) / sample_tot_vol

If you export the tables without filters (all the data from the project in 1 unique table)

Export format	Options
<input checked="" type="radio"/> General export (configurable export for general purposes)	<input checked="" type="checkbox"/> Object Data (median,mean, x, y, ...) <input checked="" type="checkbox"/> Process Data (software,version, ...) <input checked="" type="checkbox"/> Acquisition Data (Resolution, ...) <input checked="" type="checkbox"/> Sample Data (lat,long, date, ...) <input type="checkbox"/> Historical Data <input type="checkbox"/> Comments <input type="checkbox"/> Use comma as decimal separator <input type="checkbox"/> Internal Ids (including taxonomic source Id) Split in multiple files by <input type="text" value="NOT Active"/>

- You have to group the objects belonging to the same sample ->
acq_id

You have to merge d1_XXX with d2_XXX (NB : the acq_sub_part is different between d1 and d2).

- You have to group the objects belonging to the same category ->
object_annotation_category

NB: check the validation status -> object_annotation_status (have to be validated)